



The current status of the Lithium beam diagnostic at ASDEX Upgrade

Elisabeth Wolfrum
Josef Schweinzer
Matthias Reich
Rainer Fischer

Max Planck Institut für Plasmaphysik
Garching, Germany



Outline



Experimental setup

Edge ion temperatures

CX with He^{2+}

CX with D^+

Edge ion densities

Edge electron densities

ELM resolved profiles

LID evaluation with Bayes integrated concept

Wide Filters: high temporal resolution

Beam emission spectroscopy

Electron density measurement LID

Li⁰ (2p-2s) @ 670.8 nm


Li-beam: 30 – 80 keV, 2-4 mA
 \varnothing 12 mm

Charge exchange spectroscopy

Ion temperature measurement LIT

M.Reich - Thesis

Ion density measurement LIS



IslandDraw-file: /u/mat/open/diagnostics/DIV2/lib
 ../lib.eps




EM CCD improves availability of T_i measurements.





- Signal/noise ratio much improved due to EM technology
- Same radial position
- Measured temperatures agree, $T_i = 370$ eV
- Temporal resolution: 10 channels (LOS) on one CCD with 4 ms continuously

PI Micromax CCD camera




PhotonMax EM CCD camera



active CCD area




frame transfer





Extended serial register:
Electron multiplication via impact ionisation

'traditional' amplifier
1 MHz, 5 MHz

EM amplifier
5 MHz, 10 MHz



Edge ion temperature profiles




- Spatial resolution ~ 5 mm
- Temporal resolution not available:
Signal must be integrated over 1-2 s.
- He concentration $> 10\%$.
- L-mode o.k.
- H-mode: only for $f_{\text{ELM}} < 100$ Hz.




New: CX measurements also possible with D⁺ ions.




Raw data: spectrum



temporal evolution of integrated spectrum



net spectrum



- No ELMs or regular ELMs
- ELMs have to be cut out
- $\Delta t > 500$ ms integration time
- Fit difficult because centre is always dominated by photon statistics of passive line emission
- Inclusion of CXS_fit in progress




Edge ion densities, examples




C^{6+} He^{2+}



#18055, t = 3-5 s




#20463, t = 3.9-5 s





Beam emission spectroscopy BES



Lithium beam attenuation code:

$$\frac{dN_i(z)}{dz} = [n_e(z) \cdot a_{ij}(T(z)) + b_{ij}] \cdot N_j(z)$$

N_i relative occupation of state i ($i=2s, 2p, \dots, 4f, Li^+$)

a_{ij} rate coefficients

b_{ij} Einstein coefficients

a_{ij} :

Inelastic collisions with protons, electrons and impurities


b_{ij} :

radiative transitions

References:

Schweitzer et al: At.Data Nucl Data Tables 72 (1999) 239-273

Brandenburg et al: PPCF 31 (1999) 471-484



Li I Grotrian diagram



Electron density measurements




Measured profile + errors

Produce fit to data

This relative profile $Li_{2p}(z)$ is directly related to occupation number of $Li(2p)$.

$$\alpha Li_{2p}(z) = N_{2p}(z), \quad \alpha = \text{const.}$$



α is determined via 2 boundary conditions: $N_i(z=0) = \delta_{1i}$
 $N_1(z_{\text{end}}) = 0$

Use second equation of

$$\frac{dN_i(z)}{dz} = [n_e(z) \cdot a_{ij}(T(z)) + b_{ij}] \cdot N_j(z)$$

to get n_e .


Lithium beam attenuation code

Singularity near maximum of $Li(2p)$ profile.


Time interval extended to include Raus – scan:

LID with 1 ms temp. resolution


Shot: 12200 t1= 3.400 t2= 5.200



Shot: 12200 t1= 3.400 t2= 5.200



Binning of raw signal relative to ELM yields density profiles across ELM.



- Choose time interval with regular ELMs
- Determine for every t in LIB signal:
 - Δt to previous ELM
 - Δt to next ELM
- Add all signals with same temporal distance to ELM
- Calculate density profile

Attention:



ELMs should have about same size.

ELM shotfile must be checked carefully: no missed or additional ELMs.

Lithium beam must be very stable: no sparks.



n_e gradient in ETB region is recovered after 3 ms.





Determine electron density profile

given

Additional information:
 n_e profile monotonic

Description of profile:
13 knots using Hermite polynomials

Lithium beam attenuation code


Li I (2p-2s) profile

Chi² fit to the data * factor α to determine 13 knots and α .

Measured data of Li I (2p-2s)


profile and their likelihoods
! accurate error determination!

Measurements of other diagnostics.






Electron density profile evaluation now beyond point of singularity.




Medium core density:

- Pedestal well determined.
- Temporal resolution: 1 ms
- So far LID density evaluation stops just before turn.
- New: high certainty of profile up to $\rho_{\text{pol}} = 0.93$



High core density:

- So far LID density evaluation stops in SOL.
- New: reliable profile up to $\rho_{\text{pol}} = 0.96$



- New filters: 2 nm FWHM, before 0.5 nm
- Easier to change to different beam velocity (no tilt adjustment necessary)
- Higher transmission (85%, before 50%)
- Signal \sim factor 10 larger
- Now 1.5 – 4 sec with 20 kHz (before 5 kHz)



Allows faster data acquisition:
1.5 – 4 sec with 20 kHz (before 5 kHz)



Summary



Ion temperature profiles

- No temporal resolution ($\Delta t > 500$ ms)
- Good spatial resolution (5 mm)
- C: concentration now too low (< 0.4 %)
- He: good data if concentration > 10 %
- D: good data if plasma quiet, e.g. ohmic or L-mode

Ion density profiles

- Collisional mixing is important at the edge.

Electron density profiles (main business!)

- Excellent temporal resolution (~ 1 ms).
- Li-beam electron densities can resolve ELM, if several light profiles are binned relative to ELM.
- Integrated concept: edge pedestal densities can be determined up to $n_e^{\text{PED}} \sim 7 \cdot 10^{19} \text{ m}^{-3}$.
- New filters give more photons, more flexibility and allow faster data acquisition (20 kHz).